

Pediatric Positioning Chair

Optimal Positioning Yields Best Possible Images

	AP	Lateral	Oblique
• C-Spine	x	x	x
• Chest	x	x	
• Abdominal	x		
• Magnified Airway	x	x	
• Standing Lower Extremities	x		

X-rays Quickly, Accurately, and Safely

- Immobilization is done gently but securely
- Eliminates virtually all retakes, saving time and money
- The patient is secured only once – the chair rotates, positioning for different views
- Keeps patient's midline centered when moving from AP to lateral to oblique views
- Keeps patient as close to the film as possible to ensure high quality images



Fits Patients from Infants to 4 Year Olds

The uniquely designed saddle seat, with detachable saddle horn is form-fitting and very safe. The bench seat for older children includes a detachable saddle horn that aids in positioning. In either configuration, the saddle horn supports a protective genital lead shield.

- Adjustable back and head-tilt support fits patients of all sizes
- Critical chair components are x-ray transparent for excellent x-ray contrast
- Easy to keep clean, the chair is made from germ resistant, nonporous polyethylene
- High-quality lockable casters provide mobility and stability
- Simple, reliable, and virtually indestructible

Price: \$2990

Clear Image Devices' series of positioning devices is the result of innovative research and design, responding to the expressed needs of radiographers. it is our goal to produce quality products that enhance medical care and promote the continued advancement of medical imaging.